SPINTROL CAL

Multiparametric Sera / Suero Multiparamétrico Bovine Calibrator / Calibrador bovino

Serum calibrator for clinical chemistry assays IVD

Store at 2 - 8°C.

PRODUCT CHARACTERISTICS

SPINTROL CAL is a lyophilised calibrator serum from bovine source. The values given for routine methods, referred to hereafter as specified values have been selected to ensure optimum calibration of automatic analysers. Contains no additives from human source.

REAGENTS

Bovine serum. Biological additives. Bacteriostatics agents.
The concentration / activities of the components are lot-specific. The exact values and ranges valid for reagents are given informational purpose in the value sheet.

PRECAUTIONS

This calibrator is manufactured from bovine serum. Human additives which included prostatic acid phosphatase and prostatic antigen have been tested and found negative for the presence of HBsAg, HCV and antibody to HIV (1/2). However, handle cautiously as potentially infectious.

CALIBRATION

Information available on request.

PREPARATION

-Open a bottle very carefully, avoiding any loss of the lyophilized material.
-Reconstitute (→) with 5.0 mL of redistilled water.
-Mix thoroughly, avoiding foam forming.
-Bring to room temperature for about 30 min. before use.
Improper handling and/or storage can affect results.
Inaccurate reconstitution and errors in assay technique can cause erroneous results.

STORAGE AND STABILITY

The calibrator is stable until the expiration date on the label when stored tightly closed at 2-8°C and contaminations are prevented during their use. Do not use reagents over the expiration date or if there is visible evidence of microbial growth. Store calibrators tightly capped when not use.

After reconstitution is stable for:

- 15°C - 25°C 24h
- 2°C - 8°C 7 d
- 25°C - -15°C 1 month
- Bilirubin (stored protected from light):
At 2°C to 8°C 8 hours

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REAGENTS

Bovine serum. Biological additives. Bacteriostatics agents.

The concentration / activities of the components are lot-specific. The exact values and ranges valid for reagents are given informational purpose in the value sheet.

PRECAUTIONS

This calibrator is manufactured from bovine serum.

Human additives which included prostatic acid phosphatase and prostatic antigen have been tested and found negative for the presence of HBsAg, HCV and antibody to HIV (1/2). However, handle cautiously as potentially infectious.

CALIBRATION

Information available on request.

PREPARATION

-Open a bottle very carefully, avoiding any loss of the lyophilized material.

-Reconstitute (→) with 5.0 mL of redistilled water.

-Mix thoroughly, avoiding foam forming.

-Bring to room temperature for about 30 min. before use.

Improper handling and/or storage can affect results.

Inaccurate reconstitution and errors in assay technique can cause erroneous results.

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2°C - 8°C 7 d

-25°C - -15°C 1 month

- Bilirubina (stored protected from light):

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