
Human Control serum

Suero Control humano

Serum control multicomponent for clinical chemistry assays IVD

Store at 2 - 8°C.

PRODUCT CHARACTERISTICS

SPINROL H NORMAL or PATHOLOGIC is a human lyophilised serum. With most constituent concentrations and activities in the normal or pathologic range. It's intended for control of accuracy for use with manual and automated analytical procedures.

REAGENTS

Human serum. Biological additives. Bacteriostatics agents. The concentration/activities of the components are lot-specific. The exact values and ranges valid for reagents are given informational purpose in the value sheet.

PRECAUTIONS

Components from human origin have been tested and found negative for the presence of HBsAg, HCV and antibody to HIV (1/2). However handle cautiously as potentially infectious.

PREPARATION

-Reconstitute (→) with 5,0 mL of distilled water.
-Mix thoroughly, avoiding foam forming.
-Bring to room temperature for about 30 min. before use.
Improper handling and/or storage can affect results. Inaccurate reconstitution and errors in assay technique can cause erroneous results.

STORAGE AND STABILITY

The Control serum is stable until the expiration date on the label when stored tightly closed at 2-8°C and contaminations are prevented during their use. Do not use reagents over the expiration date or if there is visible evidence of microbial grown.
Store tightly capped when not use.

-After reconstitution is stable for: -Bilirubin (stored protected from light): -Acid phosphatase:
At 15°C to 25°C 12 hours At 15°C to 25°C 8 - 4 hours (total-direct) At 15°C to 25°C 4 hours
At 2°C to 8°C 5 days At 2°C to 8°C 24 - 8 hours (total-direct) At 2°C to 8°C 1 day
At -25°C to -15°C 1 month At -25°C to -15°C 2 weeks At -25°C to -15°C 2 weeks

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