HbA<sub>1c</sub>-CAL

# HbA<sub>1c</sub> Multicalibrator / Multicalibrador

Latex turbidimetry / Turbidimetria Látex

## Quantitative determination of glycated hemoglobin (HbA<sub>1c</sub>)

### IVD

Store 2 - 8°C.

### INTENDED USE

This product is for the purpose of calibrating results in the quantitative determination of human hemoglobin A1c (HbA<sub>1c</sub>) in blood by automated immunoassay

### PRODUCT CHARACTERISTICS

The lyophilized hemoglobin A1c calibrators are a hemolysate prepared from packed human erythrocytes. Stabilizers are added to maintain hemoglobin in the reduced state for the accurate calibration of the hemoglobin A1c procedure.

### REAGENTS

HbA <sub>1c</sub> Calibrator 1	Human blood. Contains preservatives. The HbA <sub>1c</sub> concentration is indicated on the vial label.
HbA <sub>1c</sub> Calibrator 2	
HbA <sub>1c</sub> Calibrator 3	
HbA <sub>1c</sub> Calibrator 4	

### PREPARATION

Reconstitute vials with 0.5mL deionized water. Gently mix for 10 minutes, or until all material has dissolved.

**The calibrators should be treated in the same manner as patient specimens regarding the hemolysate procedure.**

### PRECAUTIONS

Components from human origin have been tested and found negative for the presence of HBsAg, HCV and antibody to HIV (1/2). However handle cautiously as potentially infectious.

Avoid contact with skin and mucous membranes.

### CALIBRATION

The setpoint value of the calibrators was obtained by calibrating with materials referenced to NGSP values using the Spinreact HbA<sub>1c</sub> Reagent Set.

### MATERIAL REQUIRED AND NOT PROVIDED

Hemoglobin A1c Reagent kit, from Spinreact Ref: 43100 or 43101  
Hemoglobin A1c Control kit, from Spinreact Ref:43106

### STORAGE AND STABILITY

1. Store calibrators at 2-8°C. Stable until expiration date if sealed tightly. PROTECT FROM LIGHT AND HEAT.
2. The reconstituted calibrators should be stored refrigerated (2-8°C) and sealed tightly. The calibrator retains its assigned value for at least 30 days at 2-8°C.

### PROCEDURE

To be used in turbidimetric assays.

The lyophilized hemoglobin A1c calibrator kit will produce a calibration curve that will be stable for at least 7 days on most analyzers.

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